



AET

Disposable needle for peripheral nerve block through E.N.S. (Electro Neuro Stimulation)

A  Cannula covered by insulating material (Teflon)
Triple shaped tip
16° sharpening

B  Luer Lock connection

C  Highly insulated connection to electrostimulator

D  Extension tube

E  Electro-conducting wire with 50 cm extension



AET

Code: 0AET

Disposable needle for peripheral nerve block through E.N.S. (Electro Neuro Stimulation)

M-72 I-501 rev. A del 03/09/2012

To order:

CODE	SIZE	PACKAGE	C. CODE
0AET21070	21Gx70 mm	20	DARK GREEN
0AET21090	21Gx90 mm	20	DARK GREEN
0AET21100	21Gx100 mm	20	DARK GREEN
0AET21120	21Gx120 mm	20	DARK GREEN
0AET21150	21Gx150 mm	20	DARK GREEN
0AET22050	22Gx50 mm	20	BLACK
0AET22070	22Gx70 mm	20	BLACK
0AET22100	22Gx100 mm	20	BLACK
0AET22120	22Gx120 mm	20	BLACK
0AET22150	22Gx150 mm	20	BLACK
0AET23050	23Gx50 mm	20	DARK BLUE
0AET23070	23Gx70 mm	20	DARK BLUE
0AET23100	23Gx100 mm	20	DARK BLUE
0AET24025	24Gx25 mm	20	VIOLET
0AET24035	24Gx35 mm	20	VIOLET
0AET24050	24Gx50 mm	20	VIOLET
0AET25025	25Gx25 mm	20	ORANGE
0AET25035	25Gx35 mm	20	ORANGE
0AET25050	25Gx50 mm	20	ORANGE

On request Biomedical is available to provide needles with different lengths, tips (Chiba tip 0AEC, Quinke tip 0AEQ) and sharpening angle (45° or 30°).

biomedical

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles

Guillotine needles for biopsy on soft tissues

Cyto-aspiration needles

Different specialistic needles and sets

Loco-regional anaesthesia needles and sets

Thermic and absorbing sheets