

IM/IMX

Disposable needle for pre-operation localization of non palpable breast disease

Code: 0IM
Code: 0IMX

A



Transparent
Luer Lock

B



Indicator of the
positioned hook

C



Atraumatic special
sharpened tip
IM Version with
single hook

D




Atraumatic special
sharpened tip
IMX Version with
double hook

E



Sliding stopper
Depth marks



CE
0426



IM/IMX

Disposable needle for pre-operation localization of non palpable breast disease

Code: 0IM
Code: 0IMX

M-721-532 rev. A del 03/09/2012

To order:

CODE	SIZE	PACKAGE	C. CODE
0IM1810	18Gx10 cm	20	PINK
0IM2007	20Gx7 cm	20	YELLOW
0IM2010	20Gx10 cm	20	YELLOW
0IM2107	21Gx7 cm	20	DARK GREEN
0IM2110	21Gx10 cm	20	DARK GREEN

CODE	SIZE	PACKAGE	C. CODE
0IMX1807	18Gx7cm	20	PINK
0IMX1810	18Gx10 cm	20	PINK
0IMX2007	20Gx7 cm	20	YELLOW
0IMX2010	20Gx10 cm	20	YELLOW



On request Biomedical is available to provide needles with different lengths.

biomedical

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles
Guillotine needles for biopsy on soft tissues
Cyto-aspiration needles
Different specialistic needles and sets
Loco-regional anaesthesia needles and sets
Thermic and absorbing sheets