

MD-KIT AEQ

Kit for continuous plexus anaesthesia

Code: MD01



A 5 ml syringe
with Luer Lock cone



B AEQ needle
Quincke tip




C Catheter,
closed tip with
side holes



D Connection
between
catheter and filter



E Filter 0,22µ
0,65 ml priming volume



F Electrode



G Coupled sheet for
sterile field



CE
0373

biomedical

MD-KIT AEQ

Kit for continuous plexus anaesthesia

Code: MD01

Kit* complete with:

- I AEQ needle
- I Catheter
- I Universal adapter
- I Luer lock syringe 5 ml.
- I Antibacterial flat filter 0,2 micron
- I Polyurethane plaster
- I Electrode
- I Coupled sheet cm 50x50

*Latex free

To order:

| CODE | SIZE | C. CODE |
|------------|--------------------------------|---------|
| MD 01 0334 | 18Gx50 mm catheter 20Gx400 mm | PINK |
| MD 01 0340 | 18Gx100 mm catheter 20Gx400 mm | PINK |

On request Biomedical is available to provide needles with different lengths, tips (Chiba tip 0AEC, Tuohy tip 0AET) and sharpening angle (45° or 30°).

biomedical

distributed by: Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

| |
|---|
| Bone biopsy needles |
| Guillotine needles for biopsy on soft tissues |
| Cyto-aspiration needles |
| Different specialistic needles and sets |
| Loco-regional anaesthesia needles and sets |
| Thermic and absorbing sheets |