



MEP

Disposable needle for peripheral nerve block
under ultra sound

A  Cannula covered
by insulating material
(Teflon)
Triple shaped tip
30° sharpening

B  Luer Lock
connection

C  Extension tube



CE
0426

MEP

Code: 0MEP

Disposable needle for peripheral nerve block
under ultra sound

M-72 I-540 rev. A del 03/09/2012

To order:

CODE	SIZE	PACKAGE	C. CODE
0MEP21070	21Gx70 mm	20	DARK GREEN
0MEP21090	21Gx90 mm	20	DARK GREEN
0MEP21100	21Gx100 mm	20	DARK GREEN
0MEP21120	21Gx120 mm	20	DARK GREEN
0MEP21150	21Gx150 mm	20	DARK GREEN
0MEP22050	22Gx50 mm	20	BLACK
0MEP22070	22Gx70 mm	20	BLACK
0MEP22100	22Gx100 mm	20	BLACK
0MEP22120	22Gx120 mm	20	BLACK
0MEP22150	22Gx150 mm	20	BLACK
0MEP23050	23Gx50 mm	20	DARK BLUE
0MEP23070	23Gx70 mm	20	DARK BLUE
0MEP23100	23Gx100 mm	20	DARK BLUE
0MEP24025	24Gx25 mm	20	VIOLET
0MEP24035	24Gx35 mm	20	VIOLET
0MEP24050	24Gx50 mm	20	VIOLET
0MEP25025	25Gx25 mm	20	ORANGE
0MEP25035	25Gx35 mm	20	ORANGE
0MEP25050	25Gx50 mm	20	ORANGE

On request Biomedical is available to provide needles with different lengths, tips and sharpening angles (45° or 16°).

biomedical

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles

Guillotine needles for biopsy on soft tissues

Cyto-aspiration needles

Different specialistic needles and sets

Loco-regional anaesthesia needles and sets

Thermic and absorbing sheets