






# AV

Verres disposable needle with retractable stylet with or without stop-cock

Code: 0AVA  
Code: 0AVB  
Code: 0AVC  
Code: 0AV( $\alpha$ )Q

- A  Transparent Luer Lock  
Spring to retract stylet
- B  AVA Version with open tip window and side hole at 180°  
16° sharpening angle tip
- C  AVB Version with closed tip three windows at 120°  
16° sharpening angle tip
- D  AVC Version with closed tip window two holes at 180° to ref.  
16° sharpening angle tip
- E  AV(A)Q Needle with stop-cock mounted on Luer Lock  
30° sharpening angle tip



AV( $\alpha$ )Q  
Version with stop-cock and different sharpening angle tip

CE  
0426

biomedical



# AV

Verres disposable needle with retractable stylet with or without stop-cock

Code: 0AVA  
Code: 0AVB  
Code: 0AVC  
Code: 0AV(α)Q

M-72 I-502 rev. A del 03/09/2012

To order:

CODE	SIZE	PACKAGE	C. CODE
0AVA1410	14Gx10 cm	20	LIGHT GREEN
0AVA1416	16Gx10 cm	20	WHITE

CODE	SIZE	PACKAGE	C. CODE
0AVB1410	14Gx10 cm	20	LIGHT GREEN
0AVB1416	16Gx10 cm	20	WHITE

CODE	SIZE	PACKAGE	C. CODE
0AVC1410	14Gx10 cm	20	LIGHT GREEN
0AVC1416	16Gx10 cm	20	WHITE

CODE	SIZE	PACKAGE	C. CODE
0AV(α)Q1410	14Gx10 cm	20	LIGHT GREEN
0AV(α)Q1416	16Gx10 cm	20	WHITE



On request Biomedical is available to provide needles with different lengths and 15G diameter.

## biomedical

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy  
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629  
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles
Guillotine needles for biopsy on soft tissues
Cyto-aspiration needles
Different specialistic needles and sets
Loco-regional anaesthesia needles and sets
Thermic and absorbing sheets