

# BD

Code: 0BD

## Disposable guillotine spring-loaded needle for biopsies on soft tissues

A



Choice between two different sampling depths

B



Atraumatic sharpening tip

C



High strength steel

D



Sliding stopper  
Internal echo-marker  
Depth marks

E



Releasable external cannula with Luer Lock connection. It allows to extract several samples with one puncture only.



CE  
0426

biomedical



# BD


Code: 0BD

Disposable guillotine spring-loaded needle  
for biopsies on soft tissues

M-72 I-510 rev. A del 03/09/2012

To order:

CODE	SIZE	PACKAGE	C. CODE
0BD 1410	14Gx10 cm	20	DARK GREEN
0BD 1415	14Gx15 cm	20	DARK GREEN
0BD 1420	14Gx20 cm	20	DARK GREEN
0BD 1610	16Gx10 cm	20	WHITE
0BD 1615	16Gx15 cm	20	WHITE
0BD 1620	16Gx20 cm	20	WHITE
0BD 1810	18Gx10 cm	20	PINK
0BD 1815	18Gx15 cm	20	PINK
0BD 1820	18Gx20 cm	20	PINK
0BD 1825	18Gx25 cm	20	PINK
0BD 2010	20Gx10 cm	20	YELLOW
0BD 2015	20Gx15 cm	20	YELLOW



On request Biomedical is available to provide needles with different lengths.

## biomedical

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy  
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629  
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles

Guillotine needles for biopsy on soft tissues

Cyto-aspiration needles

Different specialistic needles and sets

Loco-regional anaesthesia needles and sets

Thermic and absorbing sheets