

BDS

Code: 0BDS

Disposable spring-loaded guillotine needle
for biopsy on soft tissues

A



Choice between
two different
sampling depths

B



Atraumatic
sharpened tip

C



Sliding stopper
Internal echo-marker
Depth marks



CE
0426

biomedical



BDS


Code: 0BDS

Disposable spring-loaded guillotine needle
for biopsy on soft tissues

M-72 I-513 rev. A del 03/09/2012

To order:

CODE	SIZE	PACKAGE	C. CODE
0BDS 1410	14Gx10 cm	20	DARK GREEN
0BDS 1415	14Gx15 cm	20	DARK GREEN
0BDS 1420	14Gx20 cm	20	DARK GREEN
0BDS 1610	16Gx10 cm	20	WHITE
0BDS 1615	16Gx15 cm	20	WHITE
0BDS 1620	16Gx20 cm	20	WHITE
0BDS 1810	18Gx10 cm	20	PINK
0BDS 1815	18Gx15 cm	20	PINK
0BDS 1820	18Gx20 cm	20	PINK
0BDS 1825	18Gx25 cm	20	PINK
0BDS 2010	20Gx10 cm	20	YELLOW
0BDS 2015	20Gx15 cm	20	YELLOW



On request Biomedical is available to provide needles with different lengths.

biomedical

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles

Guillotine needles for biopsy on soft tissues

Cyto-aspiration needles

Different specialistic needles and sets

Loco-regional anaesthesia needles and sets

Thermic and absorbing sheets