

BDPS

Code: 0BDPS

Disposable spring-loaded guillotine needle
for coaxial biopsy on soft tissues

A



Choice between
two different
sampling depths

B



Atraumatic
sharpen tip

C



Sliding stopper
Internal echo-marker
Depth marks

D



Transparent
Luer Lock

E



IT
Disposable introducer needle
for coaxial biopsy



CE
0426

BDPS


Disposable spring-loaded guillotine needle
for coaxial biopsy on soft tissues

Code: 0BDPS

M-72 I-512 rev. A del 03/09/2012

To order:

CODE	SIZE	PACKAGE	C. CODE
0BDPS 1610	16Gx10 cm	20	WHITE
0BDPS 1615	16Gx15 cm	20	WHITE
0BDPS 1620	16Gx20 cm	20	WHITE
0BDPS 1810	18Gx10 cm	20	PINK
0BDPS 1815	18Gx15 cm	20	PINK
0BDPS 1820	18Gx20 cm	20	PINK
0BDPS 2015	20Gx15 cm	20	YELLOW
0BDPS 2020	20Gx20 cm	20	YELLOW



On request Biomedical is available to provide needles with different lengths.

biomedical

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles

Guillotine needles for biopsy on soft tissues

Cyto-aspiration needles

Different specialistic needles and sets

Loco-regional anaesthesia needles and sets

Thermic and absorbing sheets