

BM

Code: 0BM

Disposable guillotine needle
for manual biopsy on soft tissues

A



Choice between
two different
sampling depths

B



Sampling depth
Limiter's reed

C




Atraumatic
sharpened tip

D



Sliding stopper
Internal echo-marker
Depth marks



CE
0426

biomedical

BM


Disposable guillotine needle
for manual biopsy on soft tissues

Code: 0BM

M-72 I-525 rev. A del 03/09/2012

To order:

CODE	SIZE	PACKAGE	C. CODE
0BM I410	14Gx10 cm	20	DARK GREEN
0BM I415	14Gx15 cm	20	DARK GREEN
0BM I420	14Gx20 cm	20	DARK GREEN
0BM I610	16Gx10 cm	20	WHITE
0BM I615	16Gx15 cm	20	WHITE
0BM I620	16Gx20 cm	20	WHITE
0BM I810	18Gx10 cm	20	PINK
0BM I815	18Gx15 cm	20	PINK
0BM I820	18Gx20 cm	20	PINK
0BM I825	18Gx25 cm	20	PINK



On request Biomedical is available to provide needles with different lengths and diameter 20G and 21G.

biomedical

 Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles
Guillotine needles for biopsy on soft tissues
Cyto-aspiration needles
Different specialistic needles and sets
Loco-regional anaesthesia needles and sets
Thermic and absorbing sheets