

BA

Disposable liver biopsy set
(according to Menghini)

Code: 0BA

A



Piston stop

B



Bisturi blade n°11

C



Cannula with
bioptic sharpening

D



Sliding stopper
Internal echo-marker
Depth marks

E



Specillum to stop
the sample

F



Hypodermic
needle

G



Luer Lock
connection



CE
0426

biomedical

BA

Disposable liver biopsy set
(according to Menghini)

Code: 0BA

To order:

CODE	SIZE	PACKAGE	C. CODE
0BA1509	15Gx9 cm	20	BLUE-GREY
0BA1609	16Gx9 cm	20	WHITE
0BA1709	17Gx9 cm	20	RED-VIOLET
0BA1809	18Gx9 cm	20	PINK

biomedical

 Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles
Guillotine needles for biopsy on soft tissues
Cyto-aspiration needles
Different specialistic needles and sets
Loco-regional anaesthesia needles and sets
Thermic and absorbing sheets