

BT/IT

Disposable introducer needle
with stylet (BT)
without stylet (IT)

Code: 0BT
Code: 0IT

A



Transparent
Luer Lock

B



Cannula with blunt
tip and stylet
with Trocar tip
(BT only)

C



Sliding stopper
Internal echo-marker
Depth marks



CE
0426

biomedical



BT/IT

Disposable introducer needle
with stylet (BT)
without stylet (IT)

Code: 0BT
Code: 0IT

M-72 I-528 rev. A del 03/09/2012

To order:

CODE	SIZE	PACKAGE	C. CODE
0BT 1410	14Gx10 cm	20	DARK GREEN
0BT 1415	14Gx15 cm	20	DARK GREEN
0BT 1420	14Gx20 cm	20	DARK GREEN
0BT 1610	16Gx10 cm	20	WHITE
0BT 1615	16Gx15 cm	20	WHITE
0BT 1620	16Gx20 cm	20	WHITE
0BT 1810	18Gx10 cm	20	PINK
0BT 1815	18Gx15 cm	20	PINK
0BT 1820	18Gx20 cm	20	PINK

CODE	SIZE	PACKAGE	C. CODE
0IT 1620	16Gx20 cm	20	WHITE
0IT 1625	16Gx25 cm	20	WHITE
0IT 1815	18Gx15 cm	20	PINK

On request Biomedical is available to provide needles with different lengths and diameters.

biomedical

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles
Guillotine needles for biopsy on soft tissues
Cyto-aspiration needles
Different specialistic needles and sets
Loco-regional anaesthesia needles and sets
Thermic and absorbing sheets