

# BE

Disposable guillotine needle for biopsies on soft tissues to be used with reusable biopsy device (Biocut-C)

A



Connectors  
for automatic devices

B



Atraumatic  
sharpened tip

C



High strength  
steel

D



Sliding stopper  
Internal echo-marker  
Depth marks



CE  
0426



# BE


Code: 0BE

Disposable guillotine needle for biopsies on soft tissues to be used with reusable biopsy device (Biocut-C)



To order:

CODE	SIZE	PACKAGE	COLOR C.
0BE 1410	14Gx10 cm	20	DARK GREEN
0BE 1415	14Gx15 cm	20	DARK GREEN
0BE 1420	14Gx20 cm	20	DARK GREEN
0BE 1610	16Gx10 cm	20	WHITE
0BE 1615	16Gx15 cm	20	WHITE
0BE 1620	16Gx20 cm	20	WHITE
0BE 1625	16Gx25 cm	20	WHITE
0BE 1810	18Gx10 cm	20	PINK
0BE 1815	18Gx15 cm	20	PINK
0BE 1820	18Gx20 cm	20	PINK
0BE 1825	18Gx25 cm	20	PINK
0BE 2015	20Gx15 cm	20	YELLOW
0BE 2020	20Gx20 cm	20	YELLOW
0BE 2025	20Gx25 cm	20	YELLOW



On request Biomedical is available to provide needles with different lengths and diameters.

## biomedical

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy  
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629  
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles

Guillotine needles for biopsy on soft tissues

Cyto-aspiration needles

Different specialistic needles and sets

Loco-regional anaesthesia needles and sets

Thermic and absorbing sheets