

# BI/BIP/BIE

Disposable needle for iliac-sternal aspiration (BI)  
Pre-assembled version (BIP)  
Explant version (BIE)

Code: 0BI  
Code: 0BIP  
Code: 0BIE




Ergonomic handle  
with fast stylet lock




Luer Lock  
connection



The threaded ring  
nut can be inserted  
or removed from  
the handle



Millimetric feed



Triple sharpened tip



BIP version  
with the pre-assembled  
grip joining the  
cannula



BIE version  
tip with lateral  
holes for explant



CE  
0426




biomedical



# BI/BIP/BIE

Disposable needle for iliac-sternal aspiration (BI)  
Pre-assembled version (BIP)  
Explant version (BIE)

Code: OBI  
Code: OBIP  
Code: OBIE



BI version mm 10/25  
BIP version mm 10/25

mm 30/45  
mm 30/45

mm 50/65  
mm 50/65

mm 75/90

To order:

## OBI

CODE	SIZE	PACKAGE	C. CODE
OBI 1410/075	14Gx10/75 mm	20	LIGHT GREEN
OBI 1425/090	14Gx25/90 mm	20	LIGHT GREEN
OBI 1510/075	15Gx10/75 mm	20	BLUE-GREY
OBI 1525/090	15Gx25/90 mm	20	BLUE-GREY
OBI 1610/075	16Gx10/75 mm	20	WHITE
OBI 1625/090	16Gx25/90 mm	20	WHITE
OBI 1810/075	18Gx10/75 mm	20	PINK
OBI 1825/090	18Gx25/90 mm	20	PINK

## OBIP (pre-assembled)

CODE	SIZE	PACKAGE	C. CODE
OBIP 1410/050	14Gx10/50 mm	20	LIGHT GREEN
OBIP 1425/065	14Gx25/65 mm	20	LIGHT GREEN
OBIP 1510/050	15Gx10/50 mm	20	BLUE-GREY
OBIP 1525/065	15Gx25/65 mm	20	BLUE-GREY
OBIP 1610/050	16Gx10/50 mm	20	WHITE
OBIP 1625/065	16Gx25/65 mm	20	WHITE
OBIP 1810/050	18Gx10/50 mm	20	PINK
OBIP 1825/065	18Gx25/65 mm	20	PINK

## OBIE

CODE	SIZE	PACKAGE	C. CODE
OBIE 1410/075	14Gx10/75 mm	20	LIGHT GREEN
OBIE 1425/090	14Gx25/90 mm	20	LIGHT GREEN
OBIE 1510/075	15Gx10/75 mm	20	BLUE-GREY
OBIE 1525/090	15Gx25/90 mm	20	BLUE-GREY
OBIE 1610/075	16Gx10/75 mm	20	WHITE
OBIE 1625/090	16Gx25/90 mm	20	WHITE
OBIE 1810/075	18Gx10/75 mm	20	PINK
OBIE 1825/090	18Gx25/90 mm	20	PINK



On request Biomedical is available to provide needles with different lengths.

# biomedical

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy  
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629  
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles

Guillotine needles for biopsy on soft tissues

Cyto-aspiration needles

Different specialistic needles and sets

Loco-regional anaesthesia needles and sets

Thermic and absorbing sheets