


BIL /BILE

Code: 0BIL
Code: 0BILE

Disposable needle for iliac-sternal aspiration
with anatomic grip (BIL) and
Explant version (BILE)




Anatomic grip




Luer Lock
connection



The threaded ring
nut can be inserted
or removed from
the handle




Millimetric feed




Triple sharpened tip



Depth limiter



BILE version
tip with lateral
holes for explant



CE
0426

biomedical




BIL /BILE

Disposable needle for iliac-sternal aspiration with anatomic grip (BIL) and Explant version (BILE)


Code: OBIL
Code: OBILE




mm 17/32




mm 36/51



mm 55/70



mm 80/95



BILE version

To order:

OBIL

CODE	SIZE	PACKAGE	C. CODE
OBIL 1417/080	14Gx17/80 mm	20	LIGHT GREEN
OBIL 1432/095	14Gx32/95 mm	20	LIGHT GREEN
OBIL 1517/080	15Gx17/80 mm	20	BLUE-GREY
OBIL 1532/095	15Gx32/95 mm	20	BLUE-GREY
OBIL 1617/080	16Gx17/80 mm	20	WHITE
OBIL 1632/095	16Gx32/95 mm	20	WHITE
OBIL 1817/080	18Gx17/80 mm	20	PINK
OBIL 1832/095	18Gx32/95 mm	20	PINK

OBILE

CODE	SIZE	PACKAGE	C. CODE
OBILE 1417/080	14Gx17/80 mm	20	LIGHT GREEN
OBILE 1432/095	14Gx32/95 mm	20	LIGHT GREEN
OBILE 1517/080	15Gx17/80 mm	20	BLUE-GREY
OBILE 1532/095	15Gx32/95 mm	20	BLUE-GREY
OBILE 1617/080	16Gx17/80 mm	20	WHITE
OBILE 1632/090	16Gx32/95 mm	20	WHITE
OBILE 1817/080	18Gx17/80 mm	20	PINK
OBILE 1832/095	18Gx32/95 mm	20	PINK

On request Biomedical is available to provide needles with different lengths.

biomedical

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles

Guillotine needles for biopsy on soft tissues

Cyto-aspiration needles

Different specialistic needles and sets

Loco-regional anaesthesia needles and sets

Thermic and absorbing sheets