



BL

Code: 0BL

Disposable needle for bone marrow biopsy version

- A  Double cup
- B  Stylet
Luer Lock connection
- C  Extraction wire
- D  Stylet with pyramidal tip
- E  Cannula with conical tightening feature to block the frustule
- F  Safety cup to extract the frustule



CE
0426

biomedical



BL

Code: 0BL

Disposable needle for bone marrow biopsy version

M-72 I-522 rev. A del 03/09/2012

To order:

CODE	SIZE	PACKAGE	C. CODE
0BL 0810	8Gx10 cm	20	DARK GREEN
0BL 0815	8Gx15 cm	20	DARK GREEN
0BL 0910	9Gx10 cm	20	PINK
0BL 0915	9Gx15 cm	20	PINK
0BL 1110	11Gx10 cm	20	GREEN YELLOW
0BL 1115	11Gx15 cm	20	GREEN YELLOW
0BL 1310	13Gx10 cm	20	LIGHT BLUE



On request Biomedical is available to provide needles with different lengths and diameter 7G.

biomedical

 Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles
Guillotine needles for biopsy on soft tissues
Cyto-aspiration needles
Different specialistic needles and sets
Loco-regional anaesthesia needles and sets
Thermic and absorbing sheets