

BLE

Disposable needle for
bone marrow explant




A
Stylet


Luer Lock
connection



B
Stylet with
pyramidal tip



C
Cannula with
lateral holes
for explant



D
Double stylet



CE
0426

BLE

Disposable needle for
bone marrow explant

Code: OBLE

M-72 I-523 rev. A del 03/09/2012

To order:

CODE	SIZE	PACKAGE	C. CODE
OBLE 0810	8Gx10 cm	20	DARK GREEN
OBLE 0815	8Gx15 cm	20	DARK GREEN
OBLE 0910	9Gx10 cm	20	PINK
OBLE 0915	9Gx15 cm	20	PINK
OBLE 1110	11Gx10 cm	20	GREEN YELLOW
OBLE 1115	11Gx15 cm	20	GREEN YELLOW
OBLE 1310	13Gx10 cm	20	LIGHT BLUE



On request Biomedical is available to provide needles with different lengths and diameter 7G.

biomedical

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles

Guillotine needles for biopsy on soft tissues

Cyto-aspiration needles

Different specialistic needles and sets

Loco-regional anaesthesia needles and sets

Thermic and absorbing sheets