








BLT

Disposable needle for bone marrow biopsy with sample trapping device

- A  Double cup
- B  Stylet
Luer Lock connection
- C  Extraction wire with depth marks
- D  Stylet with pyramidal tip
- E  Cannula with fish mouth tip and conical tightening feature to block the frustule
- F  Lock-system device to trap the sample and avoid bone biopsy deflection technique
- G  Safety cup



CE
0426



BLT

Code: OBLT

Disposable needle for bone marrow biopsy
with sample trapping device

M-72 I-524 rev. A del 03/09/2012

To order:

CODE	SIZE	PACKAGE	C. CODE
OBLT 0810	8Gx10 cm	20	DARK GREEN
OBLT 0815	8Gx15 cm	20	DARK GREEN
OBLT 0910	9Gx10 cm	20	PINK
OBLT 0915	9Gx15 cm	20	PINK
OBLT 1110	11Gx10 cm	20	GREEN YELLOW
OBLT 1115	11Gx15 cm	20	GREEN YELLOW
OBLT 1310	13Gx10 cm	20	LIGHT BLUE



On request Biomedical is available to provide needles with different lengths and diameter 7G.

biomedical

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles

Guillotine needles for biopsy on soft tissues

Cyto-aspiration needles

Different specialistic needles and sets

Loco-regional anaesthesia needles and sets

Thermic and absorbing sheets