

# EPI-KIT

Kit for continuous epidural anaesthesia

Code: EPI-KIT

A



10 ml syringe  
with Luer Lock cone

B



Tuohy needle

C



Radiopaque catheter,  
closed tip with side  
holes

D



Connection  
between  
catheter and filter

E



Filter 0,22µ  
0,65 ml priming volume



CE  
0373

biomedical

# EPI-KIT

Kit for continuous epidural anaesthesia

Code: EPI-KIT

M-721-529 rev. A del 03/09/2012

To order:

Versioni catheter without stylet

CODE	SIZE	PACKAGE	C. CODE
EPI-KIT I6	I6Gx80 mm	20	WHITE
EPI-KIT I7	I7Gx80 mm	20	RED
EPI-KIT I8	I8Gx80 mm	20	PINK

Version catheter and LOR syringe


CODE	SIZE	PACKAGE	C. CODE
EPI-KIT I6L	I6Gx80 mm	20	WHITE
EPI-KIT I7L	I7Gx80 mm	20	RED
EPI-KIT I8L	I8Gx80 mm	20	PINK

Version catheter with stylet

CODE	SIZE	PACKAGE	C. CODE
EPI-KIT I6M	I6Gx80 mm	20	WHITE
EPI-KIT I7M	I7Gx80 mm	20	RED
EPI-KIT I8M	I8Gx80 mm	20	PINK

version catheter with stylet and LOR syringe

CODE	SIZE	PACKAGE	C. CODE
EPI-KIT I6ML	I6Gx80 mm	20	WHITE
EPI-KIT I7ML	I7Gx80 mm	20	RED
EPI-KIT I8ML	I8Gx80 mm	20	PINK



**biomedical**

distributed by: Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy  
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629  
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles

Guillotine needles for biopsy on soft tissues

Cyto-aspiration needles

Different specialistic needles and sets

Loco-regional anaesthesia needles and sets

Thermic and absorbing sheets