

ID/IDK

Galactography set and kit

Code: 0ID
Code: 0IDK

ID

A




Tube
1,25x0,75x20 mm

B



Curved or straight
cannula
Transparent
Luer Lock



IDK

C



Six dilators

D



1 cc syringe

E



Curved or straight
cannula



CE
0426

biomedical



ID/IDK

Galactography set and kit

Code: 0ID
Code: 0IDK

M-72 I-53 I rev. A del 03/09/2012

To order:

CODE	SIZE	PACKAGE	C. CODE
0ID27/D20	27Gx20 mm straight cannula	20	GREY
0ID27/C20	27Gx20 mm curved cannula	20	GREY

CODE	SIZE	PACKAGE	C. CODE
0IDK25/D32	25Gx32 mm straight cannula	20	ORANGE
0IDK25/C32	25Gx32 mm curved cannula	20	ORANGE



biomedical

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles

Guillotine needles for biopsy on soft tissues

Cyto-aspiration needles

Different specialistic needles and sets

Loco-regional anaesthesia needles and sets

Thermic and absorbing sheets