

# MDI

Code: MDI

Disposable spinal atraumatic needle  
(Pencil Point)

A



Transparent cone  
to visualize the spinal  
liquid reflux

B



Complete with  
pre-assembled  
introducer

C



Atraumatic Pencil Point



CE  
0373

biomedical

# MDI

Code: MDI

## Disposable spinal atraumatic needle (Pencil Point)

### Instructions of use:

Through the needle's hole, the injected anesthetic spreads uniformly into the spinal space, obtaining an excellent result.

To order:

| CODE     | SIZE       | PACKAGE | C. CODE |
|----------|------------|---------|---------|
| MDI2290  | 22Gx90 mm  | 25      | BLACK   |
| MDI22103 | 22Gx103 mm | 25      | BLACK   |
| MDI2490  | 24Gx90 mm  | 25      | VIOLET  |
| MDI24103 | 24Gx103 mm | 25      | VIOLET  |
| MDI2590  | 25Gx90 mm  | 25      | ORANGE  |
| MDI25103 | 25Gx103 mm | 25      | ORANGE  |
| MDI25120 | 25Gx120 mm | 25      | ORANGE  |
| MDI2690  | 26Gx90 mm  | 25      | BROWN   |
| MDI26103 | 26Gx103 mm | 25      | BROWN   |
| MDI2790  | 27Gx90 mm  | 25      | GREY    |
| MDI27103 | 27Gx103 mm | 25      | GREY    |
| MDI27120 | 27Gx120 mm | 25      | GREY    |
| MDI2990  | 29Gx90 mm  | 25      | RED     |
| MDI29103 | 29Gx103 mm | 25      | RED     |

# biomedical

distributed by: Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy  
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629  
www.biomedical-bio.net sales@biomedical-bio.net

|                                               |
|-----------------------------------------------|
| Bone biopsy needles                           |
| Guillotine needles for biopsy on soft tissues |
| Cyto-aspiration needles                       |
| Different specialistic needles and sets       |
| Loco-regional anaesthesia needles and sets    |
| Thermic and absorbing sheets                  |