

MDIV

Disposable spinal atraumatic needle
(Pencil Point) with varied cross section

A



Transparent Luer Lock cone
to visualize the spinal
liquid reflux

B



Complete with
pre-assembled
introducer

C



Variation in
cross-section

D



Atraumatic Pencil Point
without sharpened edges



CE
0373

MDIV

Code: MDIV

Disposable spinal atraumatic needle
(Pencil Point) with varied cross section

Instructions of use:

The variable diameter allows a very narrow point that can be inserted causing the least possible trauma. The base of the needle, having a greater diameter than the point, has greater rigidity which permits better manoeuvrability during the spinal introduction and guarantees great safety.

M-72 I-535 rev. A del 03/09/2012

To order:

CODE	SIZE	PACKAGE	C. CODE
MDIV22/27103	22/27Gx103 mm	25	BLACK GREY

On request Biomedical is available to provide needles with different lengths.

biomedical

distributed by: Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles

Guillotine needles for biopsy on soft tissues

Cyto-aspiration needles

Different specialistic needles and sets

Loco-regional anaesthesia needles and sets

Thermic and absorbing sheets