

MI

Disposable spinal needle
for lumbar punctures

Code: MI

A



Transparent
Luer Lock

B



Quincke tip

C



Depth marks



CE
0373

biomedical



MI

Code: MI

Disposable spinal needle
for lumbar punctures

M-72 I-542 rev. A del 03/09/2012

To order:

CODE	SIZE	PACKAGE	COLOR C.
MI1890	18Gx9 cm	25	PINK
MI2090	20Gx9 cm	25	YELLOW
MI2290	22Gx9 cm	25	BLACK
MI2590	25Gx9 cm	25	ORANGE
MI2790	27Gx9 cm	25	GREY



biomedical

distributed by: Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

- Bone biopsy needles
- Guillotine needles for biopsy on soft tissues
- Cyto-aspiration needles
- Different specialistic needles and sets
- Loco-regional anaesthesia needles and sets
- Thermic and absorbing sheets