

MEF

Franseen disposable needle
for cyto-histological biopsy


Code: 0MEF

A



Transparent
Luer Lock

B



Franseen
sharpened tip

C



Sliding stopper
Internal echo-marker
Depth marks



CE
0426

biomedical

MEF

Franseen disposable needle
for cyto-histological biopsy

Code: 0MEF

M-72 I-539 rev. A del 03/09/2012

To order:

CODE	SIZE	PACKAGE	C. CODE
0MEF 1810	18Gx10 cm	20	PINK
0MEF 1815	18Gx15 cm	20	PINK
0MEF 1820	18Gx20 cm	20	PINK
0MEF 2010	20Gx10 cm	20	YELLOW
0MEF 2015	20Gx15 cm	20	YELLOW
0MEF 2020	20Gx20 cm	20	YELLOW
0MEF 2110	21Gx10 cm	20	DARK GREEN
0MEF 2115	21Gx15 cm	20	DARK GREEN
0MEF 2120	21Gx20 cm	20	DARK GREEN
0MEF 2210	22Gx10 cm	20	BLACK
0MEF 2215	22Gx15 cm	20	BLACK
0MEF 2220	22Gx20 cm	20	BLACK
0MEF 2310	23Gx10 cm	20	DARK BLUE
0MEF 2315	23Gx15 cm	20	DARK BLUE
0MEF 2320	23Gx20 cm	20	DARK BLUE

On request Biomedical is available to provide needles with different lengths.

biomedical

 Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles
Guillotine needles for biopsy on soft tissues
Cyto-aspiration needles
Different specialistic needles and sets
Loco-regional anaesthesia needles and sets
Thermic and absorbing sheets