

PG

Code: 0PG

Disposable guillotine needle for biopsy on soft tissues to be used with reusable biopsy device (Biocut-M)

A



Standard connectors for automatic devices

B



Atraumatic sharpened tip

C



High strength steel

D



Sliding stopper
Internal echo-marker
Depth marks

E



Spacer device



CE
0426

biomedical



PG

Code: 0PG

Disposable guillotine needle for biopsy on soft tissues to be used with reusable biopsy device (Biocut-M)

Technical references:


Suitable for any standard biopsy device.

Suggested to be used with Biomedial "Biocut-M" biopsy device.

M-72 I-543 rev. A del 03/09/2012

To order:

CODE	SIZE	PACKAGE	C. CODE
0PG 1410	14Gx10 cm	20	DARK GREEN
0PG 1415	14Gx15 cm	20	DARK GREEN
0PG 1420	14Gx20 cm	20	DARK GREEN
0PG 1610	16Gx10 cm	20	WHITE
0PG 1615	16Gx15 cm	20	WHITE
0PG 1620	16Gx20 cm	20	WHITE
0PG 1625	16Gx25 cm	20	WHITE
0PG 1810	18Gx10 cm	20	PINK
0PG 1815	18Gx15 cm	20	PINK
0PG 1820	18Gx20 cm	20	PINK
0PG 1825	18Gx25 cm	20	PINK
0PG 2015	20Gx15 cm	20	YELLOW
0PG 2020	20Gx20 cm	20	YELLOW
0PG 2025	20Gx25 cm	20	YELLOW



On request Biomedial is available to provide needles with different lengths and diameters.

biomedical

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

- Bone biopsy needles
- Guillotine needles for biopsy on soft tissues
- Cyto-aspiration needles
- Different specialistic needles and sets
- Loco-regional anaesthesia needles and sets
- Thermic and absorbing sheets