


SPIN EASY

Disposable light kit for
analgesia and spinal anaesthesia

Code: SPIN-EASY
Code: SPIN-V-EASY

A



Coupled sheet
for sterile field with
adhesive short side

B



Spinal
atraumatic
needle with
introducer

C



Kelly
forceps

D



0,2 micron
antibacterial
filter



Blister
package

CE
0373

biomedical

SPIN EASY

Disposable light kit for
analgesia and spinal anaesthesia

Code: SPIN-EASY
Code: SPIN-V-EASY

Kit complete with:

- I Spinal MDI needle Pencil Point with introducer*
- I Coupled sheet 75x80 cm
- I Coupled sheet 60x60 cm with hole and adhesive side
- I Kelly forceps
- I 0,2 micron antibacterial filter

M-721-544 rev. A del 03/09/2012

To order:

Version with MDI needle

CODE	SIZE	PACKAGE	C. CODE
SE 2290	22Gx90 mm	50	BLACK
SE 22103	22Gx103 mm	50	BLACK
SE 2590	25Gx90 mm	50	ORANGE
SE 25103	25Gx103 mm	50	ORANGE
SE 2790	27Gx90 mm	50	GREY
SE 27103	27Gx103 mm	50	GREY

*Version with MDIV needle variable section

CODE	SIZE	PACKAGE	C. CODE
SVE 22/27103	22/27Gx103 mm	50	GREY/BLACK

Different combinations can be realized on request

biomedical

distributed by: Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles

Guillotine needles for biopsy on soft tissues

Cyto-aspiration needles

Different specialistic needles and sets

Loco-regional anaesthesia needles and sets

Thermic and absorbing sheets