

Janus



Are you ready to open it?


Janus



Janus is the only mask in the world for the emergency ventilation support during endoscopy. It allows intra-procedural treatment of respiratory complications without interrupting the procedure. Easy and quick to apply. Safe, efficient, comfortable for the patient.


The only full-face opening mask for intra-procedural ventilation

Janus can be positioned in a few seconds around the probe already inserted, permitting the prevention or instant treatment of respiratory complications




A better control on hypoxic patient

Janus allows to treat respiratory complications without neither removing the probe nor repeat the procedure




Reduction in expenses

By applying Janus you save costs for repetition of procedure



Enhanced security

Janus reduces acute respiratory failure risk factors during endoscopic procedures



Comfort and Innovation

- Possibility to adequately sedate even fragile patients
- Better control of respiratory complications
- Operator's comfort during long procedures
- Less bother for the patient
- Enhanced security in the outpatient setting

patented

Janus



CE
0426



Fully equipped

biomedical

Biomedical Srl Via G.B. Lulli, 43 50144 Firenze Italy
Tel +39 055 352678 Fax +39 055 331629 www.biomedical-bio.net