



Snap is a new disposable device to facilitate the dispense of drugs in tablets

A Male conical connection



B Soft bag



C Removable cup



D Female conical connection




Optional  
pliers



CE

biomedical






**Snap is a new disposable device to facilitate the dispense of drugs in tablets**

Provided with the appropriate pliers, it pulverizes pills and tablets in one touch, simplifying drug's administration to patients in intensive care with a feeding tube, to patients with problems of interference and to uncooperative patients.

Made of transparent plastic, aseptical, with an innovative design and removable parts, it makes the work of doctors and nurses faster and easier.

CODE  
SNAP01

PACKAGE  
50 pcs



**biomedical**

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy  
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629  
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles
Guillotine needles for biopsy on soft tissues
Cyto-aspiration needles
Different specialistic needles and sets
Loco-regional anaesthesia needles and sets
Thermic and absorbing sheets